



Product Datasheet

TRIAD3A Antibody (orb74563)

Description	TRIAD3A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9NWF9
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart cell...



Immunohistochemical analysis of formalin...