


Product Datasheet

Avian Influenza Hemagglutinin Antibody (orb74582)

Description	Avian Influenza Hemagglutinin Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q692M2
Expiration Date	12 months from date of receipt.



The bar chart illustrates that 1ug/ml He...