



Product Datasheet

Caspase-6 Antibody (orb74600)

Description	Caspase-6 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P55212
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF7 cell lysat...



Immunohistochemical analysis of formalin...