



Product Datasheet

APG7 Antibody (orb74659)

Description	APG7 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O95352
Expiration Date	12 months from date of receipt.



Western blot analysis of Caco-2 cell lysates.



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) tissue sections.