



Product Datasheet

TLR5 Antibody (orb74671)

Description	TLR5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O60602
Expiration Date	12 months from date of receipt.



Western blot analysis
of rat brain cell ...



Immunohistochemical
analysis of formalin...