

Product Datasheet

IgM antibody (RPE) (orb746896)

Description	Goat polyclonal antibody against IgM conjugated to RPE
Species/Host	Goat
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Human IgM
Target	IgM
Preservatives	0.1% Sodium Azide (NaN ₃)Stabilizing agent (sucrose). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Concentration	Ig concentration 0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.