

## Product Datasheet

### IgA antibody (AP) (orb748223)

|                            |                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody against IgA conjugated to AP                                          |
| <b>Species/Host</b>        | Goat                                                                                           |
| <b>Reactivity</b>          | Human                                                                                          |
| <b>Conjugation</b>         | AP                                                                                             |
| <b>Tested Applications</b> | ELISA                                                                                          |
| <b>Immunogen</b>           | Human IgA.                                                                                     |
| <b>Target</b>              | IgA                                                                                            |
| <b>Preservatives</b>       | 0.1% Sodium Azide (NaN <sub>3</sub> )50% Glycerol. TRIS buffered saline, 1mM MgCl <sub>2</sub> |
| <b>Form/Appearance</b>     | Purified IgG conjugated to Alkaline Phosphatase - liquid                                       |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                                       |
| <b>Note</b>                | For research use only                                                                          |
| <b>Isotype</b>             | Polyclonal IgG                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                |