



Product Datasheet

ORAI3 Antibody (orb74838)

Description	ORAI3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BRQ5
Expiration Date	12 months from date of receipt.



Western blot analysis of A20 cell lysate...



Immunohistochemical analysis of formalin...