

Product Datasheet

IgG antibody (RPE) (orb748601)

Description	Donkey polyclonal antibody against IgG conjugated to RPE
Species/Host	Donkey
Reactivity	Rat
Conjugation	RPE
Tested Applications	FC
Immunogen	Rat IgG.
Target	IgG
Preservatives	0.1% Sodium Azide (NaN ₃)Stabilizing agent (sucrose). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.