



Product Datasheet

TDP43 Antibody (orb74863)

Description	TDP43 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13148
Expiration Date	12 months from date of receipt.



Western blot analysis of L1210 cell lysa...



Immunohistochemical analysis of formalin...