



Product Datasheet

PTK7 Antibody (orb74869)

Description	PTK7 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13308
Expiration Date	12 months from date of receipt.



Western blot analysis of (A) human colon...



Immunohistochemical analysis of formalin...