

Product Datasheet

STREPTOLYSIN-O protein (orb748726)

Description	Recombinant STREPTOLYSIN-O protein
Reactivity	Bacteria
Tested Applications	ELISA
Target	STREPTOLYSIN-O
Preservatives	None present. TRIS buffered saline
Form/Appearance	Recombinant purified Streptolysin O - liquid
Concentration	Current, batch-specific concentration 3.42 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.