

Product Datasheet

LEUKEMIA INHIBITORY FACTOR antibody (Biotin) (orb749083)

Description	Rabbit polyclonal antibody against LEUKEMIA INHIBITORY FACTOR conjugated to Biotin
Species/Host	Rabbit
Reactivity	Canine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant dog Leukemia Inhibitory Factor
Target	LEUKEMIA INHIBITORY FACTOR
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	IgG concentration 1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.