

Product Datasheet

Guinea pig IL17A protein (orb749121)

Description	Recombinant Guinea pig INTERLEUKIN-17A protein
Reactivity	Guinea pig
Tested Applications	ELISA, WB
Target	IL-17A
Preservatives	10% Trehalose. Phosphate buffered saline
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.