

Product Datasheet

CD197 antibody (RPE) (orb749135)

Description	Rat monoclonal antibody against CD197 conjugated to RPE
Species/Host	Rat
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	CCR7 rat basophilic leukaemia (RBL)-2H3 transfectants.
Target	CD197
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.