



Product Datasheet

ZBP1 Antibody (orb74914)

Description	ZBP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9H171
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse small intestine tissue.



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) tissue sections.