

Product Datasheet

HLA ABC antibody (orb749148)

Description	Mouse monoclonal antibody to HLA ABC
Species/Host	Mouse
Reactivity	Human
Conjugation	AO
Tested Applications	FC
Immunogen	Purified human tonsil lymphocyte membranes.
Target	HLA ABC
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Amethyst Orange - liquid
Concentration	IgG concentration 0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.