

Product Datasheet

CD63 antibody (CF488) (orb749336)

Description

Mouse monoclonal antibody against CD63

Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	CF488
Tested Applications	FACS, IF
Immunogen	The smooth plasma membrane fraction of MeWo cells was used as the immunogen for this CD63 antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the CD63 antibody at 2-8°C, protected from light.
Note	For research use only
Application notes	Optimal dilution of the CD63 antibody should be determined by the researcher.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P08962
Hazard Information	This CD63 antibody is available for research use only.
Dilution Range	Flow cytometry: 5ul per test per one million cells in 0.1ml or 5ul per 100ul of whole blood, Immunofluorescence: 1:50-1:100
Expiration Date	12 months from date of receipt.