

Product Datasheet

CD63 antibody (CF488) (orb749339)

Description	Mouse monoclonal antibody against CD63 conjugated to CF488
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	CF488
Tested Applications	FACS, IF
Immunogen	Full length human CD63 was used as the immunogen for this antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the CD63 antibody at 2-8°C, protected from light.
Note	For research use only
Application notes	Optimal dilution of the CD63 antibody should be determined by the researcher.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P08962
Hazard Information	This CD63 antibody is available for research use only.
Dilution Range	Flow cytometry: 5ul per test per one million cells in 0.1ml or 5ul per 100ul of whole blood, Immunofluorescence: 1:50-1:100
Expiration Date	12 months from date of receipt.