

Product Datasheet

Estrogen Receptor beta antibody (CF488) (orb749355)

Description

Mouse monoclonal antibody against

CF488

Species/Host

Mouse

Reactivity

Human

Conjugation

CF488

Tested

FACS

Applications
Immunogen

The C-terminus fragment of recombinant human Estrogen Receptor beta protein was used as the immunogen for this antibody.

Preservatives

500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the Estrogen Receptor beta antibody at 2-8°C, protected from light.

Note

For research use only

Application notes

Optimal dilution of the Estrogen Receptor beta antibody should be determined by the researcher.

Formula

500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[Q92731](#)

Hazard
Information

This Estrogen Receptor beta antibody is available for research use only.

Dilution Range

Flow cytometry: 5ul per test per one million cells in 0.1ml or 5ul per 100ul of whole blood

Expiration Date

12 months from date of receipt.