

Product Datasheet

Ku70 + Ku80 antibody (CF488) (orb749356)

Description	Mouse monoclonal antibody against Ku70 + Ku80 conjugated to CF488
Species/Host	Mouse
Reactivity	Human
Conjugation	CF488
Tested Applications	FACS, IF
Immunogen	Nuclear extract of human HL-60 cells was used as the immunogen for this Ku70 + Ku80 antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Ku70 + Ku80 antibody at 2-8°C, protected from light.
Note	For research use only
Application notes	Optimal dilution of the Ku70 + Ku80 antibody should be determined by the researcher.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Hazard Information	This Ku70 + Ku80 antibody is available for research use only.
Dilution Range	Flow cytometry: 5ul per test per one million cells in 0.1ml or 5ul per 100ul of whole blood, Immunofluorescence: 1:50-1:100
Expiration Date	12 months from date of receipt.