

Product Datasheet

von Willebrand Factor Antibody / vWF (orb749397)

Description

von Willebrand Factor (vWF) is a

cells and endothelial cells li...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

A recombinant human protein fragment was used as the immunogen for this von Willebrand Factor antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the von Willebrand Factor antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Variations in protocols, secondaries and substrates may require the von Willebrand Factor antibody to be titrated for optimal performance.1. FFPE staining requires boiling tissue sections in pH 9 10mM Tris with 1mM EDTA for 10-20 min followed by cooling at RT for 20 minutes.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G purified von Willebrand Factor antibody

Hazard

Information

This von Willebrand Factor antibody is available for research use only.

Entrez

7450

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.