

Product Datasheet

CEA antibody (CF488) (orb749427)

Description	Mouse monoclonal antibody against CEA conjugated to CF488
Species/Host	Mouse
Reactivity	Human
Conjugation	CF488
Tested Applications	IHC-P
Immunogen	Human full-length recombinant CEA protein was used as the immunogen for this antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the CEA antibody at 2-8°C, protected from light.
Note	For research use only
Application notes	Optimal dilution of the CEA antibody should be determined by the researcher.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P06731
Hazard Information	This CEA antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.