

Product Datasheet

CD45 Cocktail (Leukocyte marker) antibody (CF488) (orb749437)

Description

Mouse monoclonal antibody against

conjugated to CF488

Species/Host

Mouse

Reactivity

Human

Conjugation

CF488

Tested

FACS

Applications

Immunogen

Human B- and T-lymphoma cells were used as the immunogen for this antibody cocktail.

Preservatives

500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the CD45 antibody at 2-8°C, protected from light.

Note

For research use only

Application notes

Optimal dilution of the CD45 antibody should be determined by the researcher.

Formula

500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa + Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P08575

Hazard

Information

This CD45 antibody is available for research use only.

Dilution Range

Flow cytometry: 5ul/test/million cells in 0.1ml

Expiration Date

12 months from date of receipt.