

Product Datasheet

p27 antibody (CF488) (orb749447)

Description	Mouse monoclonal antibody against p27 conjugated to CF488
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	CF488
Tested Applications	FACS, IF
Immunogen	Recombinant human protein was used as the immunogen for the p27Kip1 antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the p27K1 antibody at 2-8°C, protected from light.
Note	For research use only
Application notes	Optimal dilution of the p27K1 antibody should be determined by the researcher.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P46527
Hazard Information	This p27K1 antibody is available for research use only.
Dilution Range	Flow cytometry: 5ul per test per one million cells in 0.1ml or 5ul per 100ul of whole blood, Immunofluorescence: 1-2ug/ml
Expiration Date	12 months from date of receipt.