



Product Datasheet

LGI2 Antibody (orb74961)

Description	LGI2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8K4Z0
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...



Immunohistochemical analysis of formalin...