

Product Datasheet

CD11a Antibody / LFA-1 (orb749611)

Description

Recognizes a protein of 180 kDa, identified as

CD11a compl...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FACS, Functional Studies, IF

Immunogen

Human CD11a protein was used as the immunogen for the CD11a antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the CD11a antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the CD11a antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P20701

Hazard Information

This CD11a antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Functional studies (order BSA/sodium azide-free format)

Expiration Date

12 months from date of receipt.