

Product Datasheet

GITR Antibody / TNFRSF18 (orb749934)

Description

GITR (glucocorticoid-induced TNFR-related

superfamily, also kn...

Species/Host

Rat

Reactivity

Mouse

Conjugation

Unconjugated

Tested Applications

FACS, IF

Immunogen

Mouse CD25+ CD4+ T cells were used as the immunogen for the GITR antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the GITR antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the GITR antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rat IgG2b, lambda

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID
035714
Hazard Information

This GITR antibody is available for research use only.

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 0.5-1ug/ml

Expiration Date

12 months from date of receipt.