

Product Datasheet

CDw17 Antibody (Lactosylceramide) (orb749936)

Description

CD17/CDw17 is an intermediate

higher gangliosides that loca...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, IF

Applications
Immunogen

Beta-2 Microglobulin associated proteins from a detergent lysate of human PBLs were used as the immunogen for the CDw17 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the CDw17 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the CDw17 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgM, kappa

Clonality

Monoclonal

Purity

PEG precipitation

Uniprot ID

Not known

Hazard
Information

This CDw17 antibody is available for research use only.

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 0.5-1ug/ml

Expiration Date

12 months from date of receipt.