

Product Datasheet

CDw60 Antibody (orb749945)

Description

This mAb reacts with 9-O-acetylated disialosyl

	is a ca...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The human rheumatoid synovial T cell line (ST-1) was used as the immunogen for the CDw60 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the CDw60 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the CDw60 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgM, kappa
Clonality	Monoclonal
Purity	PEG precipitation
Uniprot ID	Not known
Hazard Information	This CDw60 antibody is available for research use only.
Dilution Range	ELISA: 1-5ug/ml for coating (order BSA/sodium azide-free format)
Expiration Date	12 months from date of receipt.