

Product Datasheet

FOXP3 Antibody / Scurfin (orb750002)

Description

Recognizes a protein of 47-55kDa, which is

known, bu...

Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FACS, IF

Applications
Immunogen

Full-length His-tagged mouse FOXP3 protein was used as the immunogen for this FOXP3 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the FOXP3 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The optimal dilution of the FOXP3 antibody for each application should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q9BZS1

Hazard
Information

This FOXP3 antibody is available for research use only.

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 0.5-1ug/ml

Expiration Date

12 months from date of receipt.