

Product Datasheet

VLDLR Antibody / VLDL Receptor (orb750055)

Description

Very Low Density Lipoprotein Receptor

	LDL receptor gene fam...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A human recombinant partial protein from the C-terminal region was used as the immunogen for this VLDL Receptor antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the VLDL Receptor antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the VLDL Receptor antibody to be titrated up or down for optimal performance.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG1, lambda
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P98155
Hazard Information	This VLDL Receptor antibody is available for research use only.
Entrez	7436
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.