

Product Datasheet

FLI1 Antibody (orb750169)

Description

Friend leukemia integration 1 is a member of the

tha...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

Human partial recombinant protein (amino acids 107-216) was used as the immunogen for this FLI1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the FLI1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the FLI1 antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard

Information

This FLI1 antibody is available for research use only.

Entrez

[2313](#)

Dilution Range

Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.