

Product Datasheet

VLDLR Antibody / VLDL Receptor (orb750176)

Description

Very Low Density Lipoprotein Receptor

	LDL receptor gene fam...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A human recombinant partial protein from the C-terminal region was used as the immunogen for this VLDLR antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the VLDLR antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the VLDLR antibody to be titrated up or down for optimal performance.1. Staining of FFPE tissue is enhanced by boiling sections in 10mM Tris with 1mM EDTA, pH9 for 10-20 min followed by cooling at RT for 20 min.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P98155
Hazard Information	This VLDLR antibody is available for research use only.
Entrez	7436
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml (1)
Expiration Date	12 months from date of receipt

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558