

Product Datasheet

MUC2 Antibody / Mucin 2 (orb750229)

Description

Recognizes a single glycoprotein of 520kDa,

no cross-re...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Amino acids
KYPTTTPISTTTMTPTPTGTQTPTTT from the
human protein were used as the immunogen
for this Mucin-2 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin
(US sourced) and 0.05% sodium azide

Storage

Store the Mucin-2 antibody at 2-8°C (with
azide) or aliquot and store at -20°C or colder
(without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be
determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US
sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q02817

Hazard
Information

This Mucin-2 antibody is available for research
use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30
min at RT

Expiration Date

12 months from date of receipt.