


Product Datasheet

NUP107 Antibody (orb75027)

Description	NUP107 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P57740
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cell lysat...



Immunohistochemical analysis of formalin...