

Product Datasheet

Streptococcus pyogenes IdeS Protein (orb750347)

Description	IdeS from Streptococcus pyogenes, also known as Mac-1, is an endopeptidase with specificity for IgG....
Reactivity	Bacteria
Conjugation	Unconjugated
Target	IdeS
Form/Appearance	Liquid
Storage	-70°C
Tag	N-10×His
Note	For research use only
Protein Sequence	WP_010922160.1
Uniprot ID	Q9F1R7
NCBI	WP_010922160.1
Expiration Date	6 months from date of receipt.