

Product Datasheet

Human FGF R2 (IIIb) protein (orb750364)

Description	Human FGFR2 (IIIb), His Tag is expressed from human 293 cells (HEK293). It contains AA Pro 154 - Leu...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	FGF R2 (IIIb)
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.8 kDa. The protein migrates as 38-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_075259.4
Purity	90%
Source	Human FGFR2 (IIIb), His Tag (orb750364) is expressed from human 293 cells (HEK293). It contains AA Pro 154 - Leu 358 (Accession # P21802-3).
MW	24.8 kDa
Uniprot ID	P21802
NCBI	NP_075259.4
Expiration Date	6 months from date of receipt.