


Product Datasheet

IL-34 Antibody (orb75052)

Description	IL-34 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q6ZMJ4
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...