

Product Datasheet

TNFRSF4 antibody (orb750980)

Description

Mouse monoclonal antibody to TNFRSF4

Species/Host

Mouse

Reactivity

Human, Monkey

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Recombinant fragment of human OX40 (CD134) protein (around aa 59-205) (exact sequence is proprietary)

Target

TNFRSF4

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

200 µg/ml

Storage

Antibody with azide - store at 2 to 8°C.
Antibody without azide - store at -20 to -80°C.
Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

43kDa

Uniprot ID

P43489

Entrez

7293

Dilution Range

ELISA (For coating, order antibody without BSA);,Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.