

Product Datasheet

FN1 antibody (orb750990)

Description

Mouse monoclonal antibody to FN1

Species/Host	Mouse
Reactivity	Amphibian, Aves, Bovine, Equine, Human, Monkey, Porcine
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, IP, WB
Immunogen	Human plasma fibronectin
Target	FN1
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	200 µg/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	220kDa (monomer); 440kDa (dimer)
Uniprot ID	P02751
Entrez	2335
Dilution Range	Immunoprecipitation (1-2ug per 100-500ug of total protein (1ml of cell lysate)); ,Western Blot (1-2ug/ml);Immunofluorescence (1-2ug/ml); ,Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes),Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.