

Product Datasheet

TLR2 antibody (orb751011)

Description	Mouse monoclonal antibody to TLR2
Species/Host	Mouse
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant full-length human TLR2 protein.
Target	TLR2
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	200 µg/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	90kDa
Uniprot ID	O60603
Entrez	7097
Dilution Range	ELISA (For coating, order antibody without BSA);,Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.