

Product Datasheet

CEACAM5 antibody (orb751116)

Description

Mouse monoclonal antibody to CEACAM5

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF

Applications
Immunogen

Recombinant full-length human CEA protein

Target

CEACAM5

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

200 µg/ml

Storage

Antibody with azide - store at 2 to 8 °C.
Antibody without azide - store at -20 to -80 °C.
Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG2b

Clonality

Monoclonal

MW

80-200kDa

Uniprot ID

[P06731](#)

Entrez

[1048](#)

Dilution Range

Flow Cytometry (0.5-1ug/million cells);
,Immunofluorescence (1-2ug/ml),Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.