

Product Datasheet

TIA1 antibody (orb751210)

Description

Mouse monoclonal antibody to TIA1

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Recombinant human TIA1 fragment of 102 amino acid residues (aa279-380)
Target	TIA1
Preservatives	200ug/ml of Ab Purified from rabbit anti-serum by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.
Concentration	200 µg/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2b, kappa
Clonality	Monoclonal
MW	40kDa
Uniprot ID	P31483
Entrez	7072
Dilution Range	Flow Cytometry (0.5-1ug/million cells in 0.1ml); Immunofluorescence (1-2ug/ml); Western Blot (1-2ug/ml for 60 minutes at RT);Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.