

Product Datasheet

TNFRSF4 antibody (orb751271)

Description

Mouse monoclonal antibody to TNFRSF4

Species/Host

Mouse

Reactivity

Human, Monkey

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications
Immunogen

Recombinant fragment of human OX40 (CD134) protein (around aa 59-205) (exact sequence is proprietary)

Target

TNFRSF4

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

200 µg/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG2c, kappa

Clonality

Monoclonal

MW

43kDa

Uniprot ID

P43489

Entrez

7293

Dilution Range

ELISA (For coating, order antibody without BSA);,Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes),Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.