


## Product Datasheet

### GAPDH Antibody (orb75138)

|                            |                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------|
| <b>Description</b>         | GAPDH Antibody                                                                               |
| <b>Species/Host</b>        | Gallus                                                                                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                 |
| <b>Tested Applications</b> | ELISA, WB                                                                                    |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended).<br>Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                   |
| <b>Uniprot ID</b>          | <b>P04406</b>                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                              |



Western blot analysis of HeLa cell lysat...