

Product Datasheet

VIL1 antibody (orb751408)

Description

Rabbit monoclonal antibody to VIL1

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Recombinant fragment (around aa600-700) of human Villin protein (exact sequence is proprietary)

Target

VIL1

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

200 µg/ml

Storage

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Rabbit / IgG

Clonality

Monoclonal

MW

93kDa

Uniprot ID
P09327
Entrez
7429
Dilution Range

Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95 °C followed by cooling at RT for 20 minutes),Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.