



Product Datasheet

KLRA2 Antibody (orb75143)

Description	KLRA2 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q60660
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse spleen ti...



Immunohistochemical analysis of formalin...