

Product Datasheet

FOXP3/Scurfin Antibody (orb751456)

Description	FOXP3/Scurfin Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide: SQRPSRCSNPTPGP, from the C Terminus of the protein sequence according to NP_054728.2; NP_001107849.1.
Form/Appearance	Affinity Purified
Storage	Antibody can be aliquotted and stored frozen at -20°; C to -70°; C in a manual defrost freezer for six months without detectable loss of activity. The antibody can be stored at 2°; - 8°; C for 1 month without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Goat IgG
Uniprot ID	Q9BZS1
Entrez	50943
Dilution Range	IHC 2-15 ug/mL; WB 0.1-0.2 ug/mL; E 0.5-1.0 ug/mL
Expiration Date	12 months from date of receipt.