

Product Datasheet

CADM4 Antibody (orb751471)

Description	CADM4 Antibody
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Peptide with sequence C-KDERFQLEEFs, from the internal region of the protein sequence according to NP_660339.1.
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Goat IgG
Uniprot ID	Q8NFZ8
Entrez	199731
Expiration Date	12 months from date of receipt.